

Liquid Control Module Datasheet (Part Number OGI-A-1007)

The Liquid Control module is an optional accessory for the OGI Microbioreactor Mk2 Part number OGI-A-1006 to maintain the cultures at the desired optical densities, or in chemostat mode to provide a constant flow of fresh media into the culture flasks. This module contains 8 pumps for liquid handling: 4 input pumps to add fresh media into the culture flasks and 4 output pumps to remove the excess volume from the culture flasks.



Technical Specification

- Power Requirements:
 - o This unit is powered from OGI Microbioreactor Mk2 Part number OGI-A-1006
- The device is intended for indoor use only
- Optimal Performance at:
 - o Operating temperature: 15 to 40°C
 - o Relative humidity: 30% - 70%
 - o Operating altitude up to 2000 meters